clinell®

Drain Disinfectant

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Natalie Vaughan, Clinical Lead Infection Prevention & Control, Nottingham University Hospitals

Bronnen

1. Otter JA, Yezli S, French GL. The Role Played by Contaminated Surfaces in the Transmission of Nosocomial Pathogens. *Infect Control Hosp Epidemiol*. 2011;32(7):687-699. doi:10.1086/660363
2. Mitchell BG, Dancer SJ, Anderson M, Dehn E. Risk of organism acquisition from prior room occupants: A systematic review and meta-analysis. *J Hosp Infect*. 2015;91(3). doi:10.1016/j.jhin.2015.08.005
3. Weingarten RA, Johnson RC, Conlan S, et al. Genomic analysis of hospital plumbing reveals diverse reservoir of bacterial plasmids conferring carbapenem resistance. *MBio*. 2018. doi:10.1128/mBio.02011-17
4. Berrouane YF, McNutt L, Buschelman BJ, et al. Outbreak of Severe *Pseudomonas aeruginosa* Infections Caused by a Contaminated Drain in a Whirlpool Bath tub. *Clin Infect Dis*. 2002. doi:10.1086/317501
5. Grabowski M, Lobo JM, Gunnell B, et al. Characterizations of handwashing sink activities in a single hospital medical intensive care unit. *J Hosp Infect*. 2018. doi:10.1016/j.jhin.2018.04.025
6. Otter JA, Vickery K, Walker JT, et al. Surface-attached cells, biofilms and biocide susceptibility: Implications for hospital cleaning and disinfection. *J Hosp Infect*. 2015. doi:10.1016/j.jhin.2014.09.008
7. Aranega-Bou P, George RP, Verlander NQ, et al. Carbapenem-resistant Enterobacteriaceae dispersal from sinks is linked to drain position and drainage rates in a laboratory model system. *J Hosp Infect*. 2019. doi:10.1016/j.jhin.2018.12.007
8. Kotay S, Chai W, Guilford W, Barry K, Mathers AJ. Spread from the sink to the patient: In situ study using green fluorescent protein (GFP)-expressing *Escherichia coli* to model bacterial dispersion from hand-washing sink-trap reservoirs. *Appl Environ Microbiol*. 2017. doi:10.1128/AEM.03327-16
9. Mathers AJ, Vegesana K, German Mesner I, et al. Intensive care unit wastewater interventions to prevent transmission of multispecies *Klebsiella pneumoniae* carbapenemase-producing organisms. *Clin Infect Dis*. 2018. doi:10.1093/cid/ciy052
10. Deasy EC, Moloney EM, Boyle MA, et al. Minimizing microbial contamination risk simultaneously from multiple hospital washbasins by automated cleaning and disinfection of U-bends with electrochemically activated solutions. *J Hosp Infect*. 2018. doi:10.1016/j.jhin.2018.01.012
11. Muzslay M, Moore G, Alhussaini N, Wilson APR. ESBL-producing Gram-negative organisms in the healthcare environment as a source of genetic material for resistance in human infections. *J Hosp Infect*. 2017. doi:10.1016/j.jhin.2016.09.009
12. Buchan BW, Graham MB, Lindmair-Snell J, et al. The relevance of sink proximity to toilets on the detection of *Klebsiella pneumoniae* carbapenemase inside sink drains. *Am J Infect Control*. 2019. doi:10.1016/j.ajic.2018.06.021
13. Humphreys PN, Finan P, Rout S, et al. A systematic evaluation of a peracetic-acid-based high performance disinfectant. *J Infect Prev*. 2013;14(4):126-131. doi:10.1177/1757177413476125
14. Ledwoch K, Mszoga M, Pascoe M, Maillard J-Y. Dry surface biofilms: a new challenge for disinfection. In: *American Society for Microbiology*. San Francisco, United States; 2019
15. Ledwoch K, Robertson A, Laurant J, Norville P, Maillard JY. It's a trap! The development of a versatile drain biofilm model and its susceptibility to disinfection. *J Hosp Infect*. 2020;106(4):757-764. doi:10.1016/j.jhin.2020.08.010
16. Shepherd E, Leitch A, Curran E, Infection Prevention and Control Team NHS Lanarkshire. A quality improvement project to standardise decontamination procedures in a single NHS board in Scotland. *J Infect Prev*. Nov 2020;21(6):241-246
17. Garvey MI, Wilkinson MAC, Bradley CW, Holden KL, Holden E. Wiping out MRSA: effect of introducing a universal disinfection wipe in a large UK teaching hospital. *Antimicrob Resist Infect Control*. 2018;7:155

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JBN220708

GAMA Healthcare Ltd.,

The Maylands Building, Maylands Avenue, Hemel Hempstead, Hertfordshire HP2 7TG, UK.

T: +44 (0)20 7993 0030 E: info@gamahealthcare.com www.gamahealthcare.com

gama
healthcare