


**ASSUNTO:** 106 – SISTEMAS COM CONTROLADOR DE DÉBITO



<b>DESIGNAÇÃO:</b> SET À PERFUSION, FILTRE 15µ, RÉGULATEUR DE DÉBIT ET SITE D'INJECTION SANS AIGUILLE, L = 240cm PERFUSIESET, 15µ FILTER, DEBIETREGELAAR EN NAALDLOOS INSPUITPUNT, L = 240cm	<b>REFERÊNCIA</b> 110.1587	<b>DES. TÉCNICO</b> Nº1587 REV.06	<b>TIPO DE PRODUTO</b> DESCARTÁVEL
--	-------------------------------	---	---------------------------------------

**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC	C115	1
- 2-WEG BELUCHTINGSFILTER / - FILTRE À AIR	C117 + C120	1
- 15µ FILTER / - FILTRE 15 µ	C121	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 80cm, DEHP-VRIJ - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 80cm, SANS PHTHALATES (DEHP)	T678	1
-DEBIETREGELAAR / RÉGULATEUR DE DÉBIT	C108	1
- NORMALE LEIDING 140cm / - TUBULURE NORMALE 140cm	T015	1
- KLEINE KLEM / - PETITE PINCE	C43	1
- NAALDLOOS TWEERICHTINGS INSPUITPUNT - SITE D'INJECTION BIDIRECTIONNEL SANS AIGUILLE	C145	1
- TUBULURE SPECIAL 2.5 X 4.1mm, 20cm, SANS PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,1mm, 20cm, DEHP-VRIJ	T416	1
- DRAAIBARE LUER-LOCK CONNECTOR - CONNECTEUR LL MOBILE	C158 + C157	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	1



STERILE EO

<b>Embalagem</b>
Cód. Barras 560771601587
75 unidades / caixa

Dep. Quality

Ass.: José Cordeiro Date: 12 / 12 / 2019