


SUBJECT: PROLONGAMENTOS STANDARD



DESIGNAÇÃO: VERLENGLEIDING MET INGEBOUWDE TERUGSLAGKLEP, 50cm PROLONGATEUR AVEC VALVE ANTI- RETOUR INCORPORÉE, 50cm MAXIMUM PRESSURE 21bar	REFERÊNCIA 301.2544	DES. TÉCNICO Nº 2544 REV.00	TIPO DE PRODUTO DESCARTÁVEL
---	--	---	---------------------------------------

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- NKS TUBULURE 2 X 4,10mm, 50cm, SANS PHTHALATES (DEHP) - NKS LEIDING 2 X 4,10mm, 50cm, DEHP-VRIJ	T456	1
- ANTI-SIPHON VALVE	C404	1
- VROUWELIJKE PLUG MET PERFORATIE - BOUCHON FEMININ AVEC PERFORATION	C211	1

Embalagem

Cód. Barras **560771602544 1**

150 unidades / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019