


ASSUNTO: PROLONGAMENTOS STANDARD



DESIGNAÇÃO: PROLONGATEUR AVEC PERFORATEUR 2 VOIES ET LL DISTALE FEMELLE, 10 cm VERLENGLEIDING MET 2-WEGPERFORATOR EN DISTALE LL VROUWELIJK, 10 cm	REFERÊNCIA 301.2630	DES. TÉCNICO Nº 2630 REV.00	TIPO DE PRODUTO DESCARTÁVEL
--	--	---	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
2-WEG PERFORATOR / PERFORATEUR 2 VOIES	C23	1
2-WEG BELUCHTINGSFILTER / FILTRE Á AIR	C117 + C120	1
- TUBULURE NORMALE 10cm SANS PHTHALATES (DEHP) - NORMALE LEIDING 10cm – DEHP-VRIJ	T002	1
- KLEINE KLEM / - PETITE PINCE	C43	1
- VROUWELIJKE LUER-LOCK - LL FEMELLE	C45	1
- PLUG - BOUCHON	C46	1

Embalagem
Cód. Barras 560771602630 7
100 unidades / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019