


ASSUNTO: SISTEMAS COM Y

DESIGNAÇÃO: SET À PERFUSION, FILTRE 15µ, SITE D'INJECTION EN Y ET FILTRE PURGE, L = 225cm PERFUSIESET, FILTER 15µ, Y-INSPUITPUNT EN PURGEERFILTER, L = 225cm	REFERÊNCIA 110.1569	DES. TÉCNICO Nº 1569 REV.07	TIPO DE PRODUTO DESCARTÁVEL
--	--	---	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC	C115	1
- 2-WEG BELUCHTINGSFILTER / - FILTRE À AIR	C117 + C120	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- 15µ FILTER / - FILTRE 15 µ	C121	1
- TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE 190cm, SANS PHTHALATES (DEHP) - LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 190cm, DEHP-VRIJ	T653	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1
- WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- BIJSPUITPUNT IN Y, LATEXVRIJ, MET PRIKBEVEILIGING - SITE EN Y SANS LATEX, AVEC PROTECTION DE DOIGTS	C81	1
- TUBULURE SPECIAL 2.5X4.01mm, 30cm, S/PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,01mm, 30cm, DEHP-VRIJ	T419	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	1



STERILE EO

Embalagem

Cód. Barras 560771601569 1

50 unidades por saco,; 2x50= 100 unidades / caixa

Dep. Quality

Ass.: José Cordeiro Date: 12 / 12 / 2019