


ASSUNTO: TORNEIRAS E RAMPAS COM PROLONGAMENTO



DESIGNAÇÃO: RAMPE 5 ROBINETS AVEC PROLONGATEUR, 150cm KRANENBALK MET 5 KRANEN EN VERLENGLEIDING, 150cm	REFERÊNCIA 311.1516	DES. TÉCNICO Nº 1516 REV.03	TIPO DE PRODUTO DESCARTÁVEL
---	-------------------------------	--	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- KRANENBALK MET 5 KRANEN (PSU) - RAMPE 5 ROBINETS À 3 VOIES (PSU)	C134	1
- PLUG / – BOUCHON	C46	6
NKS TUBULURE 2 X 4,1mm, 150cm, SANS PHTHALATES (DEHP) NKS LEIDING 2 X 4,1mm, 150cm, DEHP-VRIJ	T459	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	2
- LUER-LOCK PLUG / – LUER-LOCK STOPJE	C28	1

Embalagem

Cód. Barras 560771601516 9

50 unidades / caixa



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019