



ASSUNTO: SISTEMAS DE SANGUE COM LUER

DESIGNAÇÃO: SET À TRANSFUSION AVEC FILTRE 200µm, 25cm², TUBULURE SANS PHTHALATES ET FILTRE PURGE TRANSFUSIESET MET FILTER 200µm, 25cm², FTALAATVRIJE LEIDING EN PURGEER FILTER	REFERÊNCIA: 206.0272 SAS25LF	DES. TÉCNICO: Nº 272 REV.03	TIPO DE PRODUTO: DESCARTÁVEL
---	---	--	--

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR, ROOD - PROTECTEUR DE PERFORATEUR ROUGE	C34	1
- PERFORATOR - PERFORATEUR	C15	1
- BLOEDTRANSFUSIEKAMER SAS 25 - CHAMBRE À TRANFUSION SAS25	C192	1
- BLOEDFILTER 25cm2, 200um - FILTRE À SANG 25cm2, 200um	C16	1
- LEIDING VOOR TRANSFUSIEKAMER - TUBULURE POUR CHAMBRE À TRANSFUSION	T206	1
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 180cm, DEHP-VRIJ - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 180cm, SANS PHTHALATES (DEHP)	T624	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1
- WIJLTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	1

Embalagem

Cód. Barras 560771600272 1

75 unidades / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019