


ASSUNTO: SISTEMAS DE SOROS ALTO DÉBITO

DESIGNAÇÃO: TROUSSE DE LAVAGE VESICAL AVEC 1 PERFORATEUR, CHAMBRE COMPTE GOUTTES ET 2 CLAMPS SET VOOR BLAASSPOELING MET 1 SPIKE, DRUPPELKAMER EN 2 KLEMMEN	REFERÊNCIA 107.2598	DES. TÉCNICO Nº 2598 REV.01	TIPO DE PRODUTO DESCARTÁVEL
---	--	---	---------------------------------------

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- 1-WEG PERFORATOR / - PERFORATEUR À UNE VOIE	C20	1
- TUBULURE 5mmx6.80mm, 25cm, SANS PHTHALATES (DEHP) - LEIDING 5mmx6.80mm, 25cm –DEHP VRIJ	T237	1
- CLAMP ON / – OFF KLEM	C19	2
- TUBULURE 5.50x7.5mm, 4.2cm, SANS PHTHALATES (DEHP) - LEIDING 5.50x7.5mm, 4.2cm –DEHP VRIJ	T240	1
- CHAMBRE, SANS PHTHALATES (DEHP) - DRUPPELKAMER DEHP VRIJ	T120	1
- TUBULURE 5.50x7.5mm-3.7cm, SANS PHTHALATES (DEHP) - LEIDING 5.50x7.5mm-3.7cm –DEHP VRIJ	T130	1
- TUBULURE 5mmx6.80mm, 170cm, SANS PHTHALATES (DEHP) - LEIDING 5mmx6.80mm, 170cm –DEHP VRIJ	T050	1
- ROLKLEM 6mm / - ROLLER 6mm	C49	1
- WIELTJE VOOR ROLKLEM 6mm / - ROULETTE POUR ROLLER	C7	1
- CONNECTOR SONDE VESICAL / BLAASSONDE	C47	1
- LEIDING 4.5x9.5mm, 10cm / - TUBULURE 4.5x9.5mm, 10cm	T591	1
- DUBBELE VERPAKKING / EMBALLAGE DOUBLE	S01	1



STERILE EO

Embalagem
Cód. Barras 560771602598 4
10 unidades por saco;
4X10= 40 unidades / caixa

Dep. Quality
Ass.: José Cordeiro Date: 10 / 10 / 2019