


ASSUNTO: SISTEMAS COM Y



| | | | |
|---|--------------------------------------|---|---------------------------------------|
| DESIGNAÇÃO: SET À PERFUSION, FILTRE 15µ, SITE D'INJECTION EN Y ET FILTRE PURGE, L = 195cm PERFUSIESET, FILTER 15µ, Y-INSPUITPUNT EN PURGEERFILTER, L = 195cm | REFERÊNCIA 110.1624 | DES. TÉCNICO Nº 1624 REV.06 | TIPO DE PRODUTO DESCARTÁVEL |
|---|--------------------------------------|---|---------------------------------------|

| COMPONENTES | | |
|---|-------------|------------|
| DESIGNAÇÃO | REFERÊNCIA | QUANTIDADE |
| - BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR | C1 | 1 |
| - PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC | C115 | 1 |
| - 2-WEG BELUCHTINGSFILTER / - FILTRE À AIR | C117 + C120 | 1 |
| - DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC | C74 | 1 |
| - 15µ FILTER / - FILTRE 15 µ | C121 | 1 |
| - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE 160cm, SANS PHTHALATES (DEHP) - LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 160cm, DEHP-VRIJ | T622 | 1 |
| - PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE | C267 | 1 |
| - WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER | C214 | 1 |
| - BIJSPUITPUNT IN Y, LATEXVRIJ, MET PRIKBEVEILIGING - SITE EN Y SANS LATEX, AVEC PROTECTION DE DOIGTS | C81 | 1 |
| - TUBULURE SPECIAL 2.5X4.01mm, 30cm, S/PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,01mm, 30cm, DEHP-VRIJ | T419 | 1 |
| - CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm | C198 | 1 |
| - LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE | C261 | 1 |



STERILE EO

Embalagem

Cód. Barras 560771601624 7

50 unidades por saco, 2x50= 100 unidades / caixa

Dep. Quality

Ass.: José Cordeiro Date: 12 / 12 / 2019