


ASSUNTO: PROLONGAMENTOS VENOSOS



DESIGNAÇÃO: VERLENGLEIDING 1x2.5mm VOOR LICHTGEVOELIGE STOFFEN, 150 cm, LL M/F PROLONGATEUR 1x2.5mm POUR PRODUITS PHOTSENSIBLES, 150cm, LL M/F	REFERÊNCIA 311.1326	DES. TÉCNICO Nº 1326 REV.01	TIPO DE PRODUTO DESCARTÁVEL
---	--------------------------------------	---	---------------------------------------

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- CONNECTOR MANNELIJK LL VOOR LEIDING 2.5mm - CONNECTEUR MALE LL POUR TUBULURE 2.5mm	C262	1
- CONNECTOR VROUWELIJK LL VOOR LEIDING 2.5mm - CONNECTEUR FEMELLE LL POUR TUBULURE 2.5mm	C263	1
- PLUG - BOUCHON	C46	1
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- LEIDING VOOR LICHTGEVOELIGE STOFFEN 1x2.5mm, 150cm, DEPH VRIJ, GEEL - TUBULURE POUR PRODUITS PHOTSENSIBLES 1x2.5mm, 150cm, SANS PHTHALATES (DEPH), JAUNE	T556	1

Embalagem

Cód. Barras 560771601326 0

50 unidades por saco,

4x50=200 unidades / caixa



STERILE EO

Dep.Quality

Ass.: José Cordeiro Date: 10 / 10 / 2019