

ASSUNTO: PROLONGAMENTOS STANDARD

DESIGNAÇÃO: PROLONGATEUR PVC, 40cm PVC VERLENGLEIDING, 40cm	REFERÊNCIA 311.0346 P2LLSP40 MF	DES. TÉCNICO Nº 346 REV.04	TIPO DE PRODUTO DESCARTÁVEL
--	--	---	---------------------------------------

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- LUER-LOCK PLUG / - LUER-LOCK STOPJE	C28	1
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- VROUWELIJKE PLUG / - BOUCHON	C46	1
- TUBULURE SPECIAL 2.5 X 4.1mm, 40cm, SANS PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,1mm, 40cm, DEHP-VRIJ	T420	1

Embalagem

Cód. Barras 560771600346 9

50 unidades por saco,

4 x 50 = 200 unidades / caixa

Maximum pressure: 24bar (350PSI)




NO DEHP



NO LATEX



STERILE EO



T420

Dep. Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019