

ASSUNTO: PROLONGAMENTOS STANDARD

DESIGNAÇÃO:	REFERÊNCIA	DES. TÉCNICO	TIPO DE PRODUTO
PROLONGATEUR PVC, 60cm PVC VERLENGLEIDING, 60cm 2.5x4.1mm	311.0348	Nº 348 REV.04	DESCARTÁVEL
	P2LLSP60 MF		

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- BESCHERMKAPJE /- BOUCHON PROTECTEUR	C28	1
- LEIDING 2.5x4.1mm, 60cm, DEPH-FREE - TUBULURE 2.5x4.1mm, 60cm, SANS PHTHALATES	T422	1
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- PLUG / - BOUCHON	C46	1

Embalagem

Cód. Barras 560771600348 3

50 unidades por saco,

4 x 50 = 200 unidades / caixa

Maximum pressure: 24bar (600PSI)



NO DEHP




NO LATEX



STERILE

EO



Dep.Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019