


ASSUNTO: SISTEMAS COM Y

DESIGNAÇÃO: SET A PERFUSION, FILTRE HYDROPHILE 8um, SITE EN Y, LUER LOCK FIXE ET FILTRE PURGE, 195cm PERFUSIESET, HYDROFIELE FILTER 8um, Y INSPUITPUNT, LL EN PURGEERFILTER, 195cm «STOPAR»	REFERÊNCIA 110.1500	DES. TÉCNICO Nº 1500 REV.06	TIPO DE PRODUTO DESCARTÁVEL
---	--	---	--



COMPONENTES

DESIGNATION	REFERENCE	RAW MATERIAL	QUANTITY
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	PE	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC	C115	ABS	1
- PRISE D'AIR AVEC FILTRE ANTIBACTÉRIEN - LUCHTINLAAT MET ANTIBACTERIELE FILTER	C290 + C120	PP / PVDF homopolymer	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	PP	1
- FILTRE HYDROPHILE 8um / - 8um HYDROFIELE FILTER	C282	Polyethersulfone	1
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 160cm, DEHP-VRIJ - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 160cm, SANS PHTHALATES (DEHP)	T622	PVC (plasticizer TOTM)	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	ABS	1
- WIJLTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	ABS	1
- BIJSPUITPUNT IN Y, LATEXVRIJ, MET PRIKBEVEILIGING - SITE EN Y SANS LATEX	C81	ABS + SEBS	1
- SPECIALE LEIDING 2.5 X 4,01mm, 30cm, DEHP-VRIJ - TUBULURE SPECIAL 2.5X4.01mm, 30cm, S/PHTHALATES (DEHP)	T419	PVC (plasticizer TOTM)	1
- LUER-LOCK CONNECTOR / - CONNECTEUR LL	C198	MABS	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	PE + Hydrophobic and Oleophobic ePTFE	1



NO DEHP



NO LATEX



STERILE EO

Dep.Quality

Ass.: José Cordeiro Date: 21/ 12 / 2020

Embalagem

Cód. Barras 560771601500 4

2x50= 100 unidades / caixa