


SUBJECT: PROLONGAMENTOS STANDARD



DESIGNAÇÃO:	REFERÊNCIA	DES. TÉCNICO	TIPO DE PRODUTO
VERLENGLEIDING MET INGEBOUWDE TERUGSLAGKLEP PROLONGATEUR AVEC VALVE ANTI-RETOUR INCORPORÉE MAXIMUM PRESSURE 21bar	301.2434	Nº 2434 REV.00	DESCARTÁVEL

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- NKS TUBULURE 2 X 4,10mm, 20cm, SANS PHTHALATES (DEHP) - NKS LEIDING 2 X 4,10mm, 20cm, DEHP-VRIJ	T489	1
- ANTI-SIPHON VALVE	C404	1
- VROUWELIJKE PLUG MET PERFORATIE - BOUCHON FEMININ AVEC PERFORATION	C211	1

Embalagem

Cód. Barras 560771602434 1

150 unidades / caixa



NO DEHP



NO LATEX



STERILE

EO

Dep.Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019