


ASSUNTO: SISTEMAS DE SOROS ALTO DÉBITO



<b>DESIGNAÇÃO:</b> TROUSSE DE LAVAGE VESICAL AVEC DOUBLE PERFORATEUR, CHAMBRE COMPTE GOUTTES ET 3 CLAMPS SET VOOR BLAASSPOELING MET 2 SPIKES, DRUPPELKAMER EN 3 KLEMMEN	<b>REFERÊNCIA</b>  <b>107.2599</b>	<b>DES. TÉCNICO</b> Nº 2599  <b>REV.01</b>	<b>TIPO DE PRODUTO</b>  DESCARTÁVEL
--	--	---	---

**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	2
- 1-WEG PERFORATOR / - PERFORATEUR Á UNE VOIE	C20	2
- TUBULURE 5mmx6.80mm, 20cm, SANS PHTHALATES (DEHP) - LEIDING 5mmx6.80mm, 20cm –DEHP VRIJ	T045	2
- CLAMP ON / – OFF KLEM	C19	3
- Y VERBINDINGSSTUK / - CONNECTION Y	C122	1
- TUBULURE 5.50x7.5mm, 4.2cm, SANS PHTHALATES (DEHP) - LEIDING 5.50x7.5mm, 4.2cm –DEHP VRIJ	T240	1
- CHAMBRE, SANS PHTHALATES (DEHP) - DRUPPELKAMER DEHP VRIJ	T120	1
- TUBULURE 5.50x7.5mm-3.7cm, SANS PHTHALATES (DEHP) - LEIDING 5.50x7.5mm-3.7cm –DEHP VRIJ	T130	1
- ROLKLEM 6mm / - ROLLER 6mm	C49	1
- WIELTJE VOOR ROLKLEM 6mm / - ROULETTE POUR ROLLER	C7	1
- CONNECTOR SONDE VESICAL / BLAASSONDE	C47	1
- LEIDING 4.5x9.5mm, 10cm / - TUBULURE 4.5x9.5mm, 10cm	T591	1
- TUBULURE 5mmx6.80mm, 170cm, SANS PHTHALATES (DEHP) - LEIDING 5mmx6.80mm, 170cm –DEHP VRIJ	T050	1
- DUBBELE VERPAKKING / EMBALLAGE DOUBLE	S01	1



STERILE EO

<b>Embalagem</b>
<b>Cód. Barras 560771602599 1</b>
30 unid. por saco, / 30 unid./ caixa

Dep.Quality  
Ass.: José Cordeiro Date: 10 / 10 / 2019