


**ASSUNTO: SISTEMAS COM TORNEIRA E Y**

<b>DESIGNAÇÃO:</b> SET À PERFUSION, FILTRE 15µ, SITE D'INJECTION EN Y,ROBINET À 3 VOIES ET FILTRE PURGE, L = 195cm PERFUSIESET, FILTER 15µ, Y-INSPUITPUNT, 3-WEGKRAAN EN PURGEERFILTER, L = 195cm	<b>REFERÊNCIA:</b>  <b>110.1610</b>	<b>DES. TÉCNICO</b> Nº 1610 <b>REV.09</b>	<b>TIPO DE PRODUTO</b> DESCARTÁVEL
---	---	---	---------------------------------------

**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- PERFORATOR PVC-VRIJ / – PERFORATEUR NON PVC	C115	1
- 2-WEG BELUCHTINGSFILTER / – FILTRE À AIR	C117 + C120	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- 15µ FILTER / – FILTRE 15 µ	C121	1
- TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE 125cm, SANS PHTHALATES (DEHP) - LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 125cm. DEHP-VRIJ	T623	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1
- WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- BIJSPUITPUNT IN Y, LATEXVRIJ, MET PRIKBEVEILIGING - SITE EN Y SANS LATEX, AVEC PROTECTION DE DOIGTS	C81	1
- TUBULURE SPECIAL 2.5X4.1mm, 30cm, S/PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,1mm, 30cm, DEHP-VRIJ	T419	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	2
- 3-WEGKRAAN / - ROBINET 3 VOIES	C260	1
- PLUG / – BOUCHON	C46	1
- NKS TUBULURE 2 X 4,1mm, 30cm, SANS PHTHALATES (DEHP) - NKS LEIDING 2 X 4,1mm, 30cm, DEHP-VRIJ	T490	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	1



<b>Embalagem</b>
<b>Cód. Barras 560771601610 0</b>
50 unidades por saco
2 x 50 = 100 unidades / caixa



STERILE EO

Dep.Quality  
Ass.: José Cordeiro Date: 12 / 12 / 2019