


ASSUNTO: PROLONGAMENTOS VENOSOS

DESIGNAÇÃO: PROLONGATEUR 1.00x3.00mm, LL M/F, 200cm VERLENGLEIDING 1.00x3.00mm LL, M/F, 200cm	REFERÊNCIA 311.1640	DES. TÉCNICO Nº 1640 REV.01	TIPO DE PRODUTO DESCARTÁVEL
--	--	---	---



COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- MANNELIJKE LUER-LOCK CONNECTOR VOOR LEIDING 3.0mm	C276	1
TUBULURE 1.00x3.00mm, 200cm, SANS PHTHALATES (DEHP) LEIDING 1.00x3.00mm, 200cm –DEHP VRIJ	T685	1
- VROUWELIJKE CYROLITE CONNECTOR 3,0	C275	1
- VROUWELIJKE PLUG / BOUCHON	C46	1

Embalagem

Cód. Barras 560771601640 7

50 unidades por saco,

3x50=150 unidades / caixa

Maximum pressure: 20bar (300PSI)

Approved by:

Ass.: José Cordeiro Data: 12 / 08 / 2019



NO DEHP



NO LATEX



STERILE EO