


**ASSUNTO: SISTEMA DE SOROS COM LUER**

<b>DESIGNAÇÃO:</b> SET POUR IRRIGATION AVEC TUBE SILICONE, VERSION COURTE IRRIGATIESET, KORTE VERSIE	<b>REFERÊNCIA</b> <b>108.1454</b>	<b>DES. TÉCNICO</b> Nº 1454 <b>REV.08</b>	<b>TIPO DE PRODUTO</b> DESCARTÁVEL
---	--------------------------------------	---	---

**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- 2-WEG PERFORATOR - PERFORATEUR 2 VOIES	C23	1
- 2-WEG BELUCHTINGSFILTER / – FILTRE À AIR	C117 + C120	1
- NORMALE LEIDING, 26cm - TUBULURE NORMALE, 26cm	T515	1
- PRECISIE ROLKLEM	C212	1
- WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- RACCORD	C66 + C112	2
- SILICONE LEIDING 3.5x5.5mm, 9cm - TUBULURE SILICONE, 3.5x5.5mm, 9cm	C84	1
- NORMALE LEIDING, 20cm - TUBULURE NORMALE, 20cm	T004	1
- VROUWELIJKE LUER-LOCK - LL FEMELLE	C45	1
- MANNELIJKE LUER-LOCK - CONNECTOR VOOR LEIDING 3,0mm	C276	1
- LEIDING 1.00x3.00mm, 220cm –DEHP VRIJ - TUBULURE 1.00x3.00mm, 220cm, SANS PHTHALATES (DEHP)	T798	1



STERILE EO

**Embalagem**

Cód. Barras 560771601454 0

100 unidades / caixa

Dep. Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019