


ASSUNTO: SISTEMAS COM TORNEIRA



<b>DESIGNAÇÃO:</b> 2 -WEG PERFUSIESET, FTALAATVRIJ FILTER 15um, 3-WEGKRAAN, DISTALE LL MET PURGEERFILTER =215cm  SET A PERFUSION AVEC CHAMBRE SANS PVC, FILTRE 15um, ROBINET 3 VOIES, LUER LOCK ET FILTRE PURGE AUTOMATIQUE =215cm	<b>REFERÊNCIA:</b>  <b>110.1868</b>	<b>DES. TÉCNICO:</b> Nº 1868  <b>REV.04</b>	<b>TIPO DE PRODUTO:</b> DESCARTÁVEL
---	---	--	--

**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC	C115	1
- 2-WEG BELUCHTINGSFILTER / - FILTRE À AIR	C117 + C120	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- 15µ FILTER / - FILTRE 15 µ	C121	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1
- WIJLTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 180cm, DEHP-VRIJ - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 180cm, SANS PHTHALATES (DEHP)	T624	1
- LUER-LOCK CONNECTOR / - CONNECTEUR LL	C33	1
- 3-WEGKRAAN / - ROBINET 3 VOIES	C260	1
- PLUG / - BOUCHON	C46	1
- NKS LEIDING 2 X 4,01mm, 15cm, DEHP-VRIJ - NKS TUBULURE 2 X 4,01mm, 15cm, SANS PHTHALATES (DEHP)	T765	1
- LUER-LOCK CONNECTOR / - CONNECTEUR LL	C198	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C261	1

Elaborado GQ

Ass.: Vânia Rodrigues Data: 20/ 10/ 2016

Embalagem

Cód. Barras 560771601868 5

50 unidades por saco,; 2x50= 100 unidades / caixa