


ASSUNTO: SISTEMAS DE SOROS COM LUER



DESIGNAÇÃO: SET A PERFUSION SANS PERFORATEUR =170cm PERFUSIESET ZONDER PERFORATOR, MAAR MET LUER LOCK AANSLUITING = 170cm	REFERÊNCIA 108.0194 S10 LL	DES. TÉCNICO Nº 194 REV.03	TIPO DE PRODUTO DESCARTÁVEL
---	---	---	---------------------------------------

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- LUER-LOCK PLUG / – LUER-LOCK STOPJE	C28	2
- LUER-LOCK CONNECTOR / – CONNECTEUR LL	C33	2
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 10cm, DEHP-VRIJ, 3.00x4.10mm - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 10cm, SANS PHTHALATES (DEHP), 3.00x4.10mm	T648	1
- DEKSEL 1-WEG VOOR DRUPPELKAMER - COUVERCLE 1 VOIE POUR CHAMBRE	C24	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 160cm, DEHP-VRIJ, 3.00x4.10mm - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 160cm, SANS PHTHALATES (DEHP), 3.00x4.10mm	T622	1
- WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1

Embalagem

Cód. Barras 560771600194 6

50 unidades por saco,

2 x 50 = 100 unidades / caixa

Elaborado GQ

Ass.: Vânia Rodrigues Data: 10 / 03 / 2016