




ASSUNTO: PROLONGAMENTOS STANDARD

DESIGNAÇÃO:	REFERÊNCIA	DES. TÉCNICO	TIPO DE PRODUTO
PROLONGATEUR PVC, 80cm PVC VERLENGLEIDING, 80cm 2.5x4.1mm	311.0349 P2LLSP80 MF	Nº 349 REV.04	DESCARTÁVEL

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- BESCHERMKAPJE /- BOUCHON PROTECTEUR	C28	1
- LEIDING 2.5x4.1mm, 80cm, DEPH-FREE - TUBULURE 2.5x4.1mm, 80cm, SANS PHTHALATES	T424	1
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- PLUG / – BOUCHON	C46	1

Embalagem
Cód. Barras 560771600349 0
50 unidades por saco,
4 x 50 = 200 unidades / caixa

Maximum pressure: 24bar (350PSI)



NO DEHP



NO LATEX



STERILE EO

Dep.Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019