


ASSUNTO: 107 - SISTEMAS DE SOROS ALTO DÉBITO

DESIGNAÇÃO: TROSSE DE LAVAGE VESICALE/ ORGANE AVEC UN PERFORATEUR ET CHAMBRE SET VOOR BLAASSPOELING, 1 SPIKE EN LL	REFERÊNCIA 107.2034	DES. TÉCNICO Nº 2034 REV.03	TIPO DE PRODUTO DESCARTÁVEL
--	--------------------------------------	---	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- 1-WEG PERFORATOR - PERFORATEUR À UNE VOIE	C20	1
- TUBULURE 5mmx6.80mm, 25cm, SANS PHTHALATES (DEHP)	T237	1
- DUBBELE KLEP	C19	1
- TUBULURE 5.50x7.5mm, 4.2cm, SANS PHTHALATES (DEHP)	T240	1
- TUBULURE 16.5mmx18mm – 14cm, SANS PHTHALATES (DEHP) - LEIDING 16.5mmx18mm – 14cm –DEHP VRIJ	T120	1
- TUBULURE 5.50x7.5mm-3.7cm, SANS PHTHALATES (DEHP)	T130	1
- TUBULURE 5mmx6.80mm, 170cm, SANS PHTHALATES (DEHP) - LEIDING 5mmx6.80mm, 170cm –DEHP VRIJ	T050	1
- ROLKLEM 6mm / - ROLLER 6mm	C49	1
- WIJLTJE VOOR ROLKLEM 6mm - ROULETTE POUR ROLLER	C7	1
- CONNECTOR	C47	1
- LEIDING 5.0x9.0mm, 10cm - TUBULURE 5.0x9.0mm, 10cm	T591	1
- LUER LOCK	C27	1
- LUER-LOCK PLUG / – LUER-LOCK STOPJE	C28	1
- DUBBELE VERPAKKING / EMBALLAGE DOUBLE	S01	1



Embalagem
Cód. Barras 560771602034 3
10 unidades por saco;
4X10= 40 unidades / caixa



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 10 / 10 / 2019