


ASSUNTO: PROLONGAMENTOS VENOSOS

DESIGNAÇÃO: PROLONGATEUR 1.00x3.00mm, LL M/F, 50cm VERLENGLEIDING 1.00x3.00mm LL, M/F, 50cm	REFERÊNCIA 311.1635	DES. TÉCNICO Nº 1635 REV.01	TIPO DE PRODUTO DESCARTÁVEL
--	--------------------------------------	---	---------------------------------------



COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- MANNELIJKE LL CONNECTOR VOOR LEIDING 3.0mm - CONNECTEUR MOLA LL POUR TUBULURE 3.0mm	C276	1
TUBULURE 1.00x3.00mm, 50cm, SANS PHTHALATES (DEHP) LEIDING 1.00x3.00mm, 50cm –DEHP VRIJ	T680	1
- VROUWELIJKE LL CONNECTOR VOOR LEIDING 3.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 3.0mm CRYOLITE	C275	1
- PLUG / BOUCHON	C46	1

Embalagem

Cód. Barras 560771601635 3

50 unidades por saco,

150 unidades / caixa

Maximum pressure: 20bar (300PSI)



NO DEHP



NO LATEX



STERILE EO

Dep.Quality

Ass.: José Cordeiro Date: 18 / 02 / 2019