


**ASSUNTO: SISTEMAS DE SOROS COM CONTROL. DÉBITO**

<b>DESIGNAÇÃO:</b> SET À PERFUSION, FILTRE HYDROPHILE, RÉGULATEUR DE DÉBIT, L = 170cm  PERFUSIESET, HYDROFIELE FILTER, DEBIETREGELAAR, L = 170cm "STOPAR"	<b>REFERÊNCIA:</b>  <b>106.2142</b>	<b>DES. TÉCNICO:</b> Nº 2142 <b>REV.01</b>	<b>TIPO DE PRODUTO:</b> DESCARTÁVEL
---	---	--	--



**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR	C1	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC	C115	1
- PRISE D'AIR AVEC FILTRE ANTIBACTÉRIEN - LUCHTINLAAT MET ANTIBACTERIELE FILTER	C290 + C120	1
- FILTRE HYDROPHILE 8um / - 8um HYDROFIELE FILTER	C282	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC	C74	1
- NORMALE LEIDING 80cm – DEHP-VRIJ - TUBULURE NORMALE 80cm S/PHTHALATES (DEHP)	T011	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE	C267	1
- WIELTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER	C214	1
- DEBIETREGELAAR / RÉGULATEUR DE DÉBIT - TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 80cm, SANS PHTHALATES (DEHP), 2.95x4.10mm - LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 80cm, DEHP-VRIJ, 2.95x4.10mm	C108	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	T678	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE	C198	1
	C261	1

**Embalagem**

**Cód. Barras 560771602142 5**

50 unidades por saco; 25 unidades por saco

50+25=75 unidades / caixa



STERILE EO

Dep.Quality  
Ass.: José Cordeiro Date: 12 / 12 / 2019