


ASSUNTO: TORNEIRAS E RAMPAS COM PROLONGAMENTO

DESIGNAÇÃO:	REFERÊNCIA	DES. TÉCNICO	TIPO DE PRODUTO
RAMPE 3 ROBINETS AVEC PROLONGATEUR, 50cm	311.1328	Nº 1328 REV.01	DESCARTÁVEL
KRANENBALK MET 3 KRANEN EN VERLENGLEIDING, 50cm	PR3PSMF50		

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- KRANENBALK MET 3 KRANEN (PSU)	C135	1
- RAMPE ROBINETS À 3 VOIES (PSU)		
- VROUWELIJKE PLUG / – BOUCHON	C46	4
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm		
- CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- TUBULURE SPECIAL 2.5 X 4.01mm, 50cm, SANS PHTHALATES (DEHP)	T421	1
- SPECIALE LEIDING 2.5 X 4,01mm, 50cm, DEHP-VRIJ		
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm		
- CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- LUER-LOCK PLUG / – LUER-LOCK STOPJE	C28	1



Embalagem

Cód. Barras 560771601328 4

50 unidades por saco,

2x50=100 / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019