


ASSUNTO: IV SET WITH LUER



DESIGNAÇÃO: SET A PERFUSION, FILTRE HYDROPHILE 8um, LUER LOCK FIXE ET FILTRE PURGE, 180cm PERFUSIESET, HYDROFIELE FILTER 8um, LL EN PURGEERFILTER, 180cm IV SET FILTER HYDROPHILIC 8um, LUER LOCK FIXE AND PURGE FILTER, 180cm	REFERÊNCIA	DES. TÉCNICO	TIPO DE PRODUTO
	108.2597	Nº 2597 REV.00	DESCARTÁVEL

COMPONENTES

DESIGNAÇÃO	REF:	QUANT
- BESCHERMKAPJE VOOR PERFORATOR - PROTECTEUR DE PERFORATEUR /- CAP	C1	1
- PERFORATOR PVC-VRIJ / - PERFORATEUR NON PVC /- SPIKE	C115	1
- PRISE D'AIR AVEC FILTRE ANTIBACTÉRIEN - LUCHTINLAAT MET ANTIBACTERIELE FILTER / - AIR INLET	C290 + C120	1
- DRUPPELKAMER PVC-VRIJ - CHAMBRE COMPTE GOUTTES SANS PVC /- CHAMBER	C74	1
- FILTRE HYDROPHILE 8um / - HYDROFIELE FILTER 8um	C282	1
- TUBULURE POUR POMPE VOLUMÉTRIQUE COMPATIBLE, 180cm, SANS PHTHALATES (DEHP), 2.95x4.10mm - LEIDING VOOR COMPATIBELE VOLUMETRISCHE POMP 180cm, DEHP-VRIJ, 2.95x4.10mm - PUMP TUBE 2.95x4.10mm, DEHP- free, 180cm	T624	1
- PRECISIE ROLKLEM ORANJE - ROLLER DE PRÉCISION ORANGE / - ORANGE ROLLER CLAMP	C267	1
- WIJLTJE VOOR PRECISIE ROLKLEM - ROULETTE POUR ROLLER / - ROLLER	C214	1
- GRANDE PINCE COULISSANTE / - GROTE KLEM - LARGE SLIDE CLAMP	C433	1
- LUER-LOCK CONNECTOR 3 X 4,1 / - CONNECTEUR LL	C198	1
- LL-PLUG MET ONTLUCHTINGSFILTER - BOUCHON LL AVEC FILTRE PURGE / - CAP WITH PURGE FILTER	C261	1

Elaborado QD
Ass.: Miguel Matias Data: 08 / 02 / 2017

Embalagem
Cód. Barras 560771602597 7
50 unidades por saco;
2x50= 100 unidades / caixa