


**ASSUNTO: PROLONGAMENTOS VENOSOS**

<b>DESIGNAÇÃO:</b> PROLONGATEUR 1.00x3.00mm, LL M/F, 10cm VERLENGLEIDING 1.00x3.00mm LL, M/F, 10cm	<b>REFERÊNCIA</b>  <b>311.1633</b>	<b>DES. TÉCNICO</b> Nº 1633 <b>REV.01</b>	<b>TIPO DE PRODUTO</b> DESCARTÁVEL
--	--	---	---



**COMPONENTES**

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1
- MANNELIJKE LL CONNECTOR VOOR LEIDING 3.0mm - CONNECTEUR MOLA LL POUR TUBULURE 3.0mm	C276	1
TUBULURE 1.00x3.00mm, 10cm, SANS PHTHALATES (DEHP) LEIDING 1.00x3.00mm, 10cm –DEHP VRIJ	T669	1
- VROUWELIJKE LL CONNECTOR VOOR LEIDING 3.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 3.0mm CRYOLITE	C275	1
- PLUG / BOUCHON	C46	1

**Embalagem**

**Cód. Barras 560771601633 9**

50 unidades por sacco,

150 unidades / caixa

**Maximum pressure: 20bar (300PSI)**



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 18 / 02 / 2019