


ASSUNTO: TORNEIRAS E RAMPAS COM PROLONGAMENTO

DESIGNAÇÃO: RAMPE 3 ROBINETS AVEC PROLONGATEUR, 25cm KRANENBALK MET 3 KRANEN EN VERLENGLEIDING, 25cm	REFERÊNCIA 311.0434 PR3PSMF25	DES. TÉCNICO Nº 1328 REV.02	TIPO DE PRODUTO DESCARTÁVEL
---	--	--	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- KRANENBALK MET 3 KRANEN (PSU) - RAMPE ROBINETS À 3 VOIES (PSU)	C135	1
- VROUWELIJKE PLUG / – BOUCHON	C46	4
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- TUBULURE SPECIAL 2.5 X 4.01mm, 25cm, SANS PHTHALATES (DEHP) - SPECIALE LEIDING 2.5 X 4,01mm, 25cm, DEHP-VRIJ	T417	1
- CONNECTOR MANNELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR MÂLE LL POUR TUBULURE 4.0mm	C198	1
- LUER-LOCK PLUG / – LUER-LOCK STOPJE	C28	1



Embalagem

Cód. Barras 560771601328 4

50 unidades por saco,

2x50=100 / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09 / 10 / 2019