


ASSUNTO: PROLONGAMENTOS STANDARD



DESIGNAÇÃO: PROLONGATEUR 15cm, 2 VALVES NON RETOUR ET LL FEMELLE VERLENGLEIDING 15cm, 2x1- RICHTINGSKLEP EN VROUWELIJKE LUER-LOCK	REFERÊNCIA 301.2691	DES. TÉCNICO Nº 2691 REV.00	TIPO DE PRODUTO DESCARTÁVEL
---	-----------------------------------	--	---

COMPONENTES

DESIGNAÇÃO	REFERÊNCIA	QUANTIDADE
- VROUWELIJKE PLUG MET PERFORATIE - BOUCHON FEMININ AVEC PERFORATION	C211	1
- CONNECTOR VROUWELIJK LL VOOR LEIDING 4.0mm - CONNECTEUR FEMELLE LL POUR TUBULURE 4.0mm	C197	1
- TUBULURE SPECIAL 2.5X4.1mm, 10cm, SANS PHTHALATES (DEHP) SPECIALE LEIDING 2.5 X 4,1mm, 10cm, DEHP-VRIJ	T414	1
- 1-RICHTINGSKLEP - VALVE NON RETOUR	C232	1
- TUBULURE SPECIAL 2.5X4.1mm, 5cm, SANS PHTHALATES (DEHP) SPECIALE LEIDING 2.5 X 4,1mm, 5cm, DEHP-VRIJ	T486	1
- 1-RICHTINGSKLEP MET LUER-LOCK - VALVE NON RETOUR	C436	1
- BESCHERMKAPJE - BOUCHON PROTECTEUR	C69	1

Embalagem

Cód. Barras 560771602691 8

100 unidades / caixa



NO DEHP



NO LATEX



STERILE EO

Dep. Quality

Ass.: José Cordeiro Date: 09/ 10 / 2019